

snap®



MILK RESIDUE TESTING

BETA-LACTAM

PENICILLIN G AMOXICILLIN AMPICILLIN CEFTIOFUR CEPHAPIRIN

TETRACYCLINE

TETRACYCLINE CHLORTETRACYCLINE OXYTETRACYCLINE

AFLATOXIN M1

AFLATOXIN M1

GENTAMICIN

GENTAMICIN

RUGGED

results in a SNAP®

SIMPLE: Self-contained and portable


CONVENIENT: No measuring of reagents, results in under 10 minutes

EASY-TO-USE: 3 quick steps with little hands-on time


3 QUICK STEPS:



Heat and pour sample




Snap device




Read

RESULTS: Clear visual read or printed results with the SNAPshot™ Reader



Negative



Positive



IDEXX
LABORATORIES

snap[®]

TESTS FROM IDEXX



- Quick screening method for fast interpretation... anywhere
- Built-in control on every test
- Consistent and reliable - trusted by dairy professionals around the world
- Results in under 10 minutes

BETA-LACTAM TEST KIT



Detects penicillin G, amoxicillin, ampicillin, ceftiofur and cephalosporin residues in raw commingled bovine milk at or below the US-FDA established tolerance/safe levels. Beta-lactams are the most widely used antibiotics in the treatment of mastitis and other infections in dairy cattle. Some common treatments which contain beta-lactam drugs are ToDAY[®], ToMORROW[®], Albacillin[®], Cefa-Lak[®], Quartermaster[®], Cefa-Dri[®], Hetacin-K[®], NAXCEL[®] and many others.

TETRACYCLINE TEST KIT

Detects tetracycline, chlortetracycline and oxytetracycline residues in raw commingled bovine milk at or below the US-FDA established tolerance/safe levels. Tetracycline antibiotics are used in the treatment of mastitis and other infections in dairy cattle. Some of the most common tetracycline treatments include LA-2000[™], Oxytet, Maxim-200[®] and many others.

AFLATOXIN M1 TEST KIT

For the screening detection of Aflatoxin M1 residues in raw commingled bovine milk.

GENTAMICIN TEST KIT

For the screening detection of gentamicin residues in raw commingled bovine milk.

SNAP[®] test kits are backed by a dedicated team of technical service professionals available to answer your product, regulatory and industry questions. Call IDEXX Technical Service at 1-800-321-0207.